

SiPM Characterization


Jeremie Cote (101069414)

Outline

- Basic theory
- Experiment goal
- Methods and results
 - Breaking voltage
 - Gain measurements
 - Spectroscopy
- Conclusion

What are SiPMs

- Silicon Photomultipliers
- Extremely sensitive photon detectors
 - Composed of a large amount of APDs
- Used more and more in high energy gamma spectroscopy and cosmic particle detection
- Array of cells can be combined to produce large detection areas




Characterization


- Breaking & bias voltage
- Gain of the SiPM
- Spectroscopy

Bias Voltage Measurements

- Measure the voltage output from the SiPM for bias voltages set on the power supply for 51 to 62 volts
- Linear fit will provide the breaking voltage and the bias voltage used will be breaking +3V



Breaking Voltage Plot



$$V_{br} = 52.268 \pm 0.699$$


$$V_p = 55.268 \pm 0.699$$

$$T\text{-Test} = 1.09$$

Gain measurement


- Photon pulses were emitted and measured by the SiPM to different channels
- Each is a plotted peak allowing for gain calculation
- Calculated by multiplying constant by distance between peaks over charge
- The overall spectrum can be Poisson fitted
- The gain for each individual set of peaks can be fit linearly against its channel
- This was done twice with the gain set to 32dB and 38dB respectively

Gain Plot




Gain 32dB = $2.634e6$

Gain 38dB = $5.353e6$



Mean +/- std dev 32dB = 7.897 +/- 3.888
Mean +/- std dev 38dB = 6.513 +/- 3.708

Gain Plot




R^2 32dB = 0.882

R^2 38dB = 0.922


Spectroscopy

- Scintillators and NaI crystal detectors used with SiPMs to measure annihilation peaks and calibrate the channels to energy


Cs-137 Spectrum



Na-22 Spectrum



Calibration



Cs-137 Energy = 0.665 ± 0.2138

Conclusion

- Overall went well, got values in the range of what was expected
- Spectroscopy for the crystal detector went very well
- Organic visibly had much poorer resolution as expected
 - Would need another sample to be able to do any estimation